

VERIFICATION OF COMPLIANCE

No.: LVD SHES2501000042HS
Applicant: MEDICAL CARRIER, SL
Calle Oficio, 20 – Pol. Ind. Los Olivos, 28906 – Getafe, Spain
Manufacturer: Same as applicant
Product Name: Pharmacy Refrigerator
Product Description: Pharmacy Refrigerator
Model No.: TH-063, TH-199
Trade Mark: 
Rating: See Page 2
Protection against Electric Shock: Class I
Additional Information: None
Sufficient samples of the product have been tested and found to be in conformity with
Test Standard: EN 61010-1:2010+A1:2019
EN IEC 61010-2-011:2021+A11:2021
as shown in the
Test Report Number(s): SHES250100004201

This Verification of Compliance has been granted to the applicant based on the results of tests, performed by Laboratory of SGS-CSTC Standards Technical Services Co., Ltd. on sample of the above-mentioned product in accordance with the provisions of the relevant harmonized standards under the Low Voltage Directive 2014/35/EU. The CE marking can be affixed, under the responsibility of the manufacturer, after completion of an EC Declaration of Conformity and compliance with all relevant EC Directives. The affixing of the CE marking presumes in addition that the conditions in annexes III and IV of the Directive are fulfilled.

Andrew Zhao

Technical Manager

SGS-CSTC Standards Technical Services Co., Ltd.

2025-01-08

This document is issued by the Company under its General Conditions of Service accessible at http://www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

SGS-CSTC Standards Technical Services Co., Ltd.
16F, Century Yu Hui Mansion, No.73, Fucheng Road, Haidian District, Beijing, 100142, CHINA
www.sgsgroup.com.cn
Member of the SGS Group (SGS SA)

No.:

LVD SHES2501000042HS

Other information added:

Rating: 220 V - 240 VAC, 50 Hz;

Model	Refrigerant	Rated current	PS(High/Low)
TH-063	R600a / 16 g	0,9 A	0,98 Mpa / 0,3 Mpa
TH-199	R600a / 25 g	1,8 A	0,65 Mpa / 0,1 Mpa

Andrew Zhao

Technical Manager

SGS-CSTC Standards Technical Services Co., Ltd.

2025-01-08

This document is issued by the Company under its General Conditions of Service accessible at http://www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

SGS-CSTC Standards Technical Services Co., Ltd.
16F, Century Yu Hui Mansion, No.73, Fucheng Road, Haidian District, Beijing, 100142, CHINA
www.sgsgroup.com.cn
Member of the SGS Group (SGS SA)